pastelle PRO

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

660nm


3mm

Possible to **remove multiple-color of tattoos** using dye

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



Principles of treatment using Genesis

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[色素沈着と組合せ] Maximize **pigmented lesion removal effect** with **chromophore pre-heating**

[トーニングと組合せ] Maximize **collagen production & skin revitalization effects** and **toning effect**

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




Low Recurrence Rate

* PTP: Photoacoustic Twin Pulse

Kim, J. Y et al. (2016). Treatment of melasma with the photoacoustic twin pulse mode of low-fluence 1,064 nm Q-switched Nd: YAG laser. *Annals of dermatology*, 28(3), 290-296.

Yoo, M. G. et al.(2015). Influence of pulse type on subcellular selective photothermolysis of melanosomes in adult zebrafish skin following 1,064-nm, Q-switched, Nd: YAG laser irradiation: a pilot study. *Annals of Dermatology*, 27(2), 230-232.

Multiple PTP vs. Conventional Toning

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Principles of treatment using Multiple PTP **pastelle PRO**

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


pISSN 1013-9087 · eISSN 2005-3894
http://dx.doi.org/10.5021/ad.2016.28.3.290

ORIGINAL ARTICLE

Treatment of Melasma with the Photoacoustic Twin Pulse Mode of Low-Fluence 1,064 nm Q-Switched Nd:YAG Laser

Jee Young Kim, Misoo Choi, Chan Hee Nam, Ji Seok Kim, Myung Hwa Kim, Byung Cheol Park, Seung Phil Hong

Department of Dermatology, Dankook University College of Medicine, Cheonan, Korea

Recently, new Q-switched Nd:YAG lasers which are capable of utilizing the photoacoustic twin pulse (PTP) mode, such as RevLite® (HOYA ConBio, Fremont, CA, USA) and Pastelle® (Won Technology Co., Ltd, Daejeon, Korea), have been developed for high-efficiency laser toning. PTP is the latest technology in which successive two beams are emitted at very short intervals (100~130 μ sec) and then synergistically produce higher peak power than those produced by a single beam from the current Nd:YAG lasers. Although there were some previous reports regarding low-fluence Nd:YAG laser therapy for melasma, little information exists about the efficacy and safety of laser toning with PTP mode. In this study, we performed an open-label controlled clinical experiment to evaluate the degree of improvement and side effects of PTP mode 1,064 nm Q-switched Nd:YAG laser toning in Asian women with moderate-to-severe melasma. Using the results from previous studies as historical controls, this study was aimed to investigate the efficacy and safety of laser toning with PTP mode.

treatment. Subjects were advised to avoid direct sun exposure and to apply broad-spectrum sunscreen during and after the treatment.

Efficacy assessment

Clinical and photographic evaluations for improvement (a decrease in the appearance of melasma) were conducted at baseline and at the last follow-up visit (2 weeks after the last treatment). High-quality digital photographs were taken from the front and sides under the same lighting conditions using a digital camera (D400; Nikon, Tokyo, Japan). Responses to treatments were evaluated at the last follow-up visit by using MASI scoring, colorimeter measurement, and the investigators' and patients' overall assessments. The MASI values were determined to quantify changes in pigmentation, as previously described^{12,13}. In addition, objective measurement of skin color was conducted at the darkest spot with a colorimeter (CR-400; Minolta, Tokyo, Japan). Measurement of brightness (L* value) were obtained three times, and the average value

DISCUSSION

Q-switched Nd:YAG lasers emit a beam at the level of nanoseconds pulse width to selectively destroy melanosomes without affecting surrounding tissue. Especially, the collimated 1,064 nm Q-switched Nd:YAG laser which delivers laser beam evenly throughout the dermis is able to treat a pigmented lesion located deep within the dermis. For achieving stable treatment results in melasma, it is important to minimize the unwanted side effects such as inflammation or skin irritation. Laser toning with 1,064 nm Q-switched Nd:YAG laser was devised as a way to reduce rebound aggravation of melasma or post-inflammatory hyperpigmentation due to such secondary stimulations which can be caused by other existing treatment modalities. The mechanism of laser toning has been proposed to be due to the fragmentation of subcellular-specific organelle melanin granules in the dendrites and subsequent dispersion into the cytoplasm without cellular destruction by repetitive doses of laser energy with a sub-photothermolytic fluence (<5 J/cm²) over a large spot size, the so called "subcellular selective photothermolysis"^{14,15}.

PTP is a unique laser emitting mode in which double pulses are delivered within one Q-switching cycle. Each pulse has relatively weak energy compared to the standard single Q-switched beam, but they can transfer higher peak energy (up to 60% more) to the target melanosome than a single beam because double beams are successively irradiated at very short intervals (100~130 μ sec) and their energy can be accumulated. As a result, the high peak energy instantly increases the temperature of the chromophore, leading to pressure changes and vibration, which then effectively destroys the chromophore through the form of a shock wave (photoacoustic effect).

Q-switching engineering in this study produces a pulse width of 20 ns on single Q-switching. The pulse width between the first and the second Q-switching in the double Q-switching PTP mode is about 100 μ s. The instant peak energy of the singly Q-switched beam is 1.3 J, whereas the energy of doubly Q-switched beams (each beam energy is about 0.8 J) in the PTP mode is 1.6 J, and as a result, the instant peak energy is increased up to 123% (<http://www.wtlaser.com/kr/products/pastelle.asp>). PTP can

months to obtain satisfactory clinical improvement^{5,11}. When compared with the results of other studies using low fluence Q-switched Nd:YAG treatments, even after considering methodological variability, about 60% of the subjects in our study had 'significantly improved', whereas previous publications showed an improvement range of 50%~90%^{2,4,5,10,16}. Therefore, the degree of improvement of global assessment in our study was not likely to be superior. However, it should be considered that five sessions of laser toning monotherapy were conducted in this study and it is a fewer treatment session than those of other previous studies, which conducted 5 to 15 times of mono- or combined-laser therapies. Moreover, the objective colorimeter L* value increased by 1.3 point in our study, and this showed better improvement as compared to an increase of 1.1 point after laser toning therapy of 10 times weekly in the same ethnic group¹⁷.

Similar to this study, a prospective single arm clinical study was recently conducted in which 50 Korean women with melasma were treated with 15 weeks of weekly treatments, using another PTP mode Q-switched Nd:YAG laser at 1,064 nm with a fluence of 2.8 J/cm², and an 8-mm spot size⁴. They reported that both patients and investigators rated the treatment outcome as "good improvement" on average with improvement rate of $\geq 50\%$ and <75% without any serious side effects. Considering the fewer treatment sessions, the more increased level of skin brightness value by using objective colorimeter after laser monotherapy, the efficacy of PTP laser toning in our study is thought to be of a significant level. It can be presumed that further increasing the number of treatment sessions may show a better effect.

Interestingly, the overall global assessment of patients revealed more improvement than that of the investigators, and this can be assumed to be due to the low incidence of discomfort or side effects and the feeling that the skin is generally clearer after the treatment. In fact, a large number of subjects were pleased with the results and noted their facial hue became brighter through the treatments (data not shown). When assessing the MASI score at the final visit, the improvement of dark color was more prominent than the degree of regression of melasma area (data not shown). As the treatment progressed, eventually, the

Irradiation by the beams dividing one pulse to two within 100-130 μ s produces higher energy than with the generally used pulse

As a result, the instant peak energy at PTP mode can be transferred to skin up to 60% more than at the standard mode, which then is more successful to destroy the chromophore

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WONTECH Co., Ltd.